

**Quantitative Determination of GLUCOSE in Serum/Plasma**  
Only for *In Vitro* Diagnostic use

**ORDER INFORMATION**

REF	Cont.
GLUSLR 500	1 X 500 ML
GLUSLR 500	5 X 100 ML
GLUSLR 1000	1 X 1000 ML

**CLINICAL SIGNIFICANCE**

Glucose is a major source of energy for most cells of the body; insulin facilitates glucose entry into the cells. Diabetes is a disease manifested by hyperglycemia patients with diabetes demonstrate an inability to produce insulin. Clinical diagnosis should not be made on a single test result; it should integrate clinical and other laboratory data.

**PRINCIPLE**

Glucose is determined after enzymatic oxidation in the presence of glucose oxidase. The hydrogen peroxide formed reacts, under catalysis of peroxidase, with phenol and 4-aminophenazone to form a red-violet quinonemine dye as indicator.

**REAGENT COMPOSITION**

Reagent I : Glucose Reagent  
Glucose standard : 100 mg/dl (5.54 mm)

**SAMPLE COLLECTION AND PRESERVATION**

Serum, plasma, and urine.  
Serum and plasma stability: 7 days at +2 -8°C

**REAGENT PREPARATION**

The Reagent is ready to use.

**REAGENT STORAGE AND STABILITY**

Reagent is stable till expiry when stored at 2 - 8°C.  
Store protected from light.

AUTOMATED PARAMETERS	
Wavelength	505 nm (490-550 nm)
Reaction Type	End Point
Cuvette	1 cm light path
Reaction Temperature	37°C
Reaction Type	Increasing
Measurement	Against Reagent Blank
Sample Volume	10µl
Reagent Volume	1000µl
Incubation	05 minutes
Blank Absorbance Limit	< 0.300
Low Normal at 37°C	60mg/dl (3.33 mmol/L)
High Normal at 37°C	110 mg/dl (6.1 mmol/L)
Linearity at 37°C	400mg/dl (22.2 mmol/L)

**MANUAL ASSAY PROCEDURE**

**PIPETTE INTO TEST TUBES**

	BLANK	STD	SAMPLE
Sample	-	-	10µl
Standard	-	10µl	-
Reagent	1000µl	1000µl	1000µl

Mix & Incubate for 05 min. at 37°C or 15 min. at R.T. Measure absorbance of Sample (AT) and Standard (AS) against Reagent Blank at 505 nm. The colour is stable for 30 min. at R.T.

**CALCULATION**

$$\text{Total Glucose (mg/dl)} = \text{AT/AS} \times \text{Conc of Standard}$$

**LINEARITY**

The method is linear upto a concentration of 400 mg/dl (22.2 mmol/L). Dilute samples above this concentration 1:1 and multiply the result by 2.

**REFERENCE VALUES**

Serum/Plasma (Fasting)	70 - 110 mg/dl
(2 hrs. P. P.)	upto 150 mg/dl
Urine	<0.5 g/24 h

**INTERFERENCES**

**Haemoglobin:** No interference up to 1000 mg/dL.  
**Free Bilirubin:** No interference up to 787 µmol/L (46 mg/dL).  
**Conjugated Bilirubin:** No interference up to 189 µmol/L (11 mg/dL).  
**Lipaemia:** No interference from lipaemia, measured as triglycerides up to 23 mmol/L (2000 mg/dL).  
**Ascorbate:** No interference from ascorbate up to 0.36 mmol/L (6.5 mg/dL)..

**QUALITY CONTROL**

It is recommended to run a normal and a pathological control serum which is commercially available to verify the performance of the measured procedure. The value of controls should fall within the established limit.

**Notes**

In case of urine, samples should be diluted 1:10 and the result multiplied by 10.

**BIBLIOGRAPHY**

Teitz, N.W., Fundamentals of Clin. Chem., Philadelphia. W.B. Saunders (1970) Trinder.P., "Determination of Blood Glucose Using 4 Aminophenazone."