

**Quantitative determination of Creatinine in serum & Urine**  
**Only for *In Vitro* Diagnostic use**

**ORDER INFORMATION**

REF	Cont.
CRZ 125	5 X 25 ML
CRZ 375	15 X 25 ML
CRZ 600	10 X 60 ML
CRZ 40	1 X 40 ML
CRZ 80	1 X 80 ML

**CLINICAL SIGNIFICANCE**

Creatinine is the catabolic product of high energy storage compound, Creatinine Phosphate formed in muscle. The amount of creatinine produced is fairly constant and is primarily a function of muscle mass. Creatinine is excreted out of body entirely by the kidneys. Elevated levels are found in renal dysfunction, reduced renal blood flow (shock, dehydration, congestive heart failure) diabetes acromegaly. Decreased levels are found in muscular dystrophy.

**Method**

Photometric Enzymatic Test method

**PRINCIPLE**

Through a series of enzymatic reactions, creatinine is converted in to glycine, whilst endogenous components such as creatine and sarcosine are eliminated in the first step of the sequence. The formed hydrogen peroxide reacts with TOPS in the presence of peroxidase, to give a quinoneimine dye. The intensity of color, measured at 546 nm, is proportional to creatinine concentration in the sample.

**REAGENT**

Reagent I : Buffer Reagent  
Reagent II : Enzyme Reagent  
Creatinine Standard : 2 mg/dl (0.16 mmol/L)

**REAGENT PREPARATION**

Use separate reagents ready to use.

**REAGENT STORAGE AND STABILITY**

Stability: up to expiration date on labels at 2-8°C.  
Stability since first opening of vials: use preferably within 60 days at 2-8°C.

**Prior to use:** When stored at 2-8°C and protected from direct sunlight, the reagents are stable until the expiry date stated on the bottle and kit box labels.

**WARNING AND PRECAUTIONS**

- For in vitro diagnostic use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Exercise the normal precautions required for handling all laboratory reagents.
- The reagent contains preservative. Do not swallow. Avoid contact with skin and mucous membranes.
- For detailed information refer Material Safety Data Sheet.

**WASTE MANAGEMENT**

Please refer to local legal requirements.

**MATERIALS REQUIRED BUT NOT PROVIDED**

- NaCl solution 9 g/L
- General laboratory equipment

**SAMPLE COLLECTION AND PRESERVATION**

**Serum or Urine**

It is very important to store the sample protected from light!  
Stability: 1 day at 20 – 25°C

7 days at 4 – 8°C  
6 months at –20°C in case of immediate freezing.  
in Urine: 1 day at 20 – 25°C  
4 days at 4 – 8°C  
3 weeks at –20°C  
Freeze only once! Discard contaminated specimens!

**ASSAY PROCEDURE**

**Operating Instructions**

- Check reagent inventories at least daily to ensure that quantities are sufficient for the planned work load.
- Bring all reagents, standard and samples to room temperature 18 – 28 °C, prior to analysis.

AUTOMATED PARAMETERS	
Wavelength	650nm (620 – 650nm)
Cuvette Light Path	1 cm
Reaction Type	End Point
Reaction Temperature	37°C
Measurement	Against Reagent blank
Sample Volume	10 µl
Reagent 1 Volume	300 µl
Reagent 2 Volume	100 µl
Incubation	5 mins. + 5 mins.
Low Normal	0.80 mg/dl
High Normal	1.40 mg/dl
Linearity	25 mg/dl

**MANUAL ASSAY PROCEDURE**

**Pipette into Test Tubes**

Dispense:	Blank	Standard	Sample
Reagent R1	300 µl	300 µl	300 µl
Water	10 µl	-	-
Calibrator	-	10 µl	-
Sample	-	-	10 µl
<b>Mix, incubate at 37°C for 5 minutes.</b>			
Dispense:	Blank	Standard	Sample
Reagent R2	100 µl	100 µl	100 µl
<b>Mix, incubate at 37°C for 5 minutes.</b>			
Read absorbance at 546 nm of Standard (AS) and samples (AT) against reagent blank.			

**SAMPLE DILUTIONS**

- This method is linear upto a concentration of 25 mg/dL.
- Dilute samples above this concentration 1:1 with 0.9% saline
- Repeat assay. Multiply the result by 2.

**CALCULATION**

**Serum/plasma sample:**

Creatinine mg/dl = AT/AS x Conc of Standard

**Random urine sample:**

Creatinine mg/dl = AT/AS x Conc of Standard x 20

**24 hours urine sample (Creatinine mg/24h):**

Creatinine mg/24h = AT/AS x Conc of Standard x 20 x diuresis (dilution, diuresis in dL)

**CALIBRATORS AND CONTROLS**

For the calibration of automated photometric systems the commercially available suitable multi-calibrator is recommended.

It is recommended to run a normal and a pathological control serum which is commercially available to verify the performance of the measured procedure. The value of controls should fall within the established limit.

Each laboratory should establish corrective action in case of deviations in control recovery.

## PERFORMANCE CHARACTERISTICS

### WITHIN RUN

Sample	Mean Concentration	SD	CV %
Norm	1.03	0.02	2.34%
Path	3.35	0.09	2.60%

### RUN TO RUN

Sample	Mean Concentration	SD	CV %
Norm	1.00	0.03	1.67%
Path	3.44	0.06	1.66%

### LINEARITY

The method is linear upto a concentration of 25mg/dL. Dilute samples above this concentration 1:1 with 0.9% saline solution and repeat assay. Multiply the result by 2.

**Limit of detection:** The limit of detection for Creatinine Total is 0.04 mg/dL.

### METHOD COMPARISON

A comparison of Accucare Bilirubin Total with a commercially available assay (x) using 20 samples gave following results:  $R^2 = 0.9900$

### REFERENCE VALUES

SPECIMEN	MEN	WOMEN
SERUM	0.8 - 1.4	0.7 - 1.2 mg/dl
24h URINE	1.0 - 2.0	0.8 - 1.8 G/24h

The reference values are to be considered as indicative only. Every laboratory should establish its own normal range.

### LIMITATION OF THE PROCEDURE

- For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.











### INTERFERENCE


- Hemoglobin: No interference found upto 1000 mg/dL.
- Lipemia: No Interference found upto 1400mg/dL.
- Bilirubin: No Interference found upto 28mg/dL.
- Ascorbic Acid: No Interference found upto 50mg/dL.
- These characteristics have been obtained using an automatic analyzer. Results may vary if a different instrument or a manual procedure is used.

### BIBLIOGRAPHY

- Tietz Textbook of Clinical Chemistry, Fourth Edition, Burtis-Ashwood-Bruns (2006), 797-801
- Clin. Chem. 2012, 58(2), 391-401

### GLOSSARY OF SYMBOL

	Consult Instruction for Use		Lot Number
	Catalog Number		Date of Manufacturing
	Store between		Use By or Expiration Date
	Manufacturer		For <i>in vitro</i> Diagnostic use only
	Keep away from sunlight		Content of the kit

 LAB-CARE DIAGNOSTICS (INDIA) PVT. LTD.  
C1 Type, Shed No.: 3225, Chemical Zone,  
GIDC Sarigam – 396155, Dist. Valsad, Gujarat, India.  
Tel.: +91 22 2554 2109 /1558  
Email: accucarediagnostics.com; Website: www.labcarediagnostics.com