

Quantitative determination of ADA activity in human Serum / Plasma / other body fluids
Only for *In Vitro* Diagnostic use

ORDER INFORMATION

REF	Pack Size
ADA 15	1 X 15 mL
ADA 30	1 X 30 mL

CLINICAL SIGNIFICANCE


ADA is an enzyme catalyzing the deamination reaction from adenosine to inosine. The enzyme is widely distributed in human tissues, and is especially high in T lymphocytes. Elevated serum ADA activity has been observed in patients with acute hepatitis, alcoholic hepatic fibrosis, chronic active hepatitis, liver cirrhosis, viral hepatitis and hepatoma. Increased ADA activity was also observed in patients with tuberculous effusions. Determination of ADA activity in patient serum may add unique values to the diagnosis of liver diseases in combination with ALT or γ -GT (GGT) tests. ADA assay may also be useful in the diagnosis of tuberculous pleuritis.

METHOD

Photometric Enzymatic Test.

PRINCIPLE

The ADA assay is based on the enzymatic deamination of adenosine to inosine which is converted to hypoxanthine by purine nucleoside phosphorylase (PNP). Hypoxanthine is then converted to uric acid and hydrogen peroxide (H₂O₂) by xanthine oxidase (XOD). H₂O₂ is further reacted with N-Ethyl-N-(2-hydroxy-3-sulfopropyl)-3-methylaniline (TOOS) and 4-aminoantipyrine (4-AA) in the presence of peroxidase (POD) to generate quinone dye which is monitored in a kinetic manner. The entire enzymatic reaction scheme is shown below.



One unit of ADA is defined as the amount of ADA that generates one μ mole of inosine from adenosine per min. at 37°C.

REAGENT

Reagent 1 : Enzyme solution
Reagent 2 : Substrate solution
ADA Calibrator : Concentration is stated on Vial label

MATERIALS REQUIRED BUT NOT PROVIDED

- NaCl solution 9 g/L
- General laboratory equipment

REAGENT PREPARATION

The reagents are provided in a ready to use format.

REAGENT STORAGE AND STABILITY

Reagents are stable until the expiry date shown on the label when stored tightly closed at 2 – 8 °C and if contaminations are prevented during their use

WARNING AND PRECAUTIONS

- For in vitro diagnostic use.
- Do not use components beyond the expiration date.
- Do not mix materials from different kit lot numbers.
- Exercise the normal precautions required for handling all laboratory reagents.
- The reagent contains preservative. Do not swallow. Avoid contact with skin and mucous membranes.

- For detailed information refer Material Safety Data Sheet.

WASTE MANAGEMENT

Please refer to local legal requirements.

MATERIALS REQUIRED BUT NOT PROVIDED

- NaCl solution 9 g/L
- General laboratory equipment

SAMPLE COLLECTION AND PRESERVATION

Serum or heparin plasma or other body fluids (Pleural Fluid, Peritoneal fluid, Pericardial fluid and CSF)

Stability:

7 days at 4 – 8 °C

1 month at –20 °C

Discard contaminated specimens. Only freeze once!

ASSAY PROCEDURE

Operating Instructions

- Check reagent inventories at least daily to ensure that quantities are sufficient for the planned work load.
- Bring all reagents, standard and samples to room temperature 18 – 28°C, prior to analysis.

Automated Parameters	
Wavelength	546 nm
Measurement	Against DI Water
Reaction Temperature	37°C
Reaction Type	Fix time Kinetic
Reaction Direction	Increasing
Incubation	3 Min.
Sample Volume	10 μ l
Reagent I Volume	400 μ l
Reagent II Volume	200 μ l
Delay/Lag/Time	300 Sec.(5min.)
Measuring Time	180 Sec.(3min.)
Linearity	200
Units	IU/L

MANUAL ASSAY PROCEDURE

Pipette into Test Tubes

REAGENT I	400 μ l
SAMPLE	10 μ l
Mix well and incubate for 3 mins at 37°C & Immediately Add	
REAGENT II	200 μ l
Mix & aspirate. After the initial delay time of 300 seconds, read the abs.of the test for the next 180 seconds at 546 nm. Determine the change in absorbance and calculate the test result.	

SAMPLE DILUTIONS

- The method is linear to a concentration of 200 IU/L.
- If the concentration exceeds this value, the sample should be diluted 1:1 with 0.9% saline solution and reassayed. Multiply the result by 2.

CALCULATION

ADA activity (IU/L) = $\frac{\Delta O.D./min. \text{ of sample}}{\Delta O.D./min. \text{ of Calibrator}} \times \text{Calibrator Conc.}$

CALIBRATORS AND CONTROLS

For the calibration of automated photometric systems the commercially available suitable multi-calibrator is recommended.

This method has been traceable to the reference materials ERM-BCR 647.

It is recommended to run a normal and a pathological control serum which is commercially available to verify the performance of the measured procedure. The value of controls should fall within the established limit.

PERFORMANCE CHARACTERISTICS WITHIN RUN

Sample	Mean Concentration	SD	CV %
Normal	31.77	1.87	5.89%
Path	90.87	2.04	2.25%

RUN TO RUN

Sample	Mean Concentration	SD	CV %
Normal	31.58	1.59	5.06%
Path	92.63	3.17	3.29%

LINEARITY

The method is linear to a concentration of 200 IU/L. If the concentration exceeds this value, the sample should be diluted 1:1 with 0.9% saline solution and reassayed. Multiply the result by 2.

Limit of detection: The limit of detection for Adenosine Deaminase is 4IU/L.

METHOD COMPARISON

A comparison of Accucare Adenosine Deaminase with a commercially available assay (x) using 20 samples gave following results: $R^2 = 0.9800$

REFERENCE VALUES

for serum, plasma, Pleural paracardial and Ascitic Fluid

Normal up to 43 U/L
Suspect for MTB 43 to 62 U/L
Strong Suspect for MTB greater than 62 U/L

For CSF

Normal up to 11 U/L
Suspect for TBM 11 to 12.35 U/L
Strong Suspect for TBM greater than 12.35 U/L
MTB (Mycobacterium Tuberculosis)
TBM (Tuberculous Meningitis)

Each laboratory should check if the reference ranges are transferable to its own patient population and determine own reference ranges if necessary.

LIMITATION OF THE PROCEDURE

- For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings.







INTERFERENCE


- Ascorbic Acid: No interference found upto 50mg/dl.
- Hemoglobin: No interference found upto Hemoglobin 800mg/dl.
- Lipemia: No interference found upto 1000mg/dl.
- These characteristics have been obtained using an automatic analyzer. Results may vary if a different instrument or a manual procedure is used.

BIBLIOGRAPHY

- Kobayashi F, Ikeda T, Marumo F, Sato C: Adenosine deaminase isoenzymes in liver disease. Am.J.Gastroenterol. 88:266-271(1993)
- Kalkan A., Bult V., Erel O., Avci S., and Bingol N.K.: Adenosine Deaminase and guanosine deaminase activities in sera of patients with viral hepatitis. Mem Inst. Oswaldo Cruz 94 (3) 383-386 (1999).

GLOSSARY OF SYMBOL

	Consult instruction for Use	LOT	Lot Number
REF	Catalog Number		Date of Manufacturing
	Store between		Use By or Expiration Date
	Manufacturer	IVD	For <i>in vitro</i> Diagnostic use only
	Keep away from sunlight	CONT	Content of the kit

 LAB-CARE DIAGNOSTICS (INDIA) PVT. LTD.
C1 Type, Shed No.: 3225, Chemical Zone,
GIDC Sarigam – 396155, Dist. Valsad, Gujarat, India.
Tel.: +91 22 2554 2109 /1558
Email: accucarediagnostics.com; Website: www.labcarediagnostics.com