

INTENDED USE

For the Quantitative determination of micro albumin (μ ALB)
Only for *In Vitro* Diagnostic use

ORDER INFORMATION

| REF | CONT |
|---------|------------|
| MAL 50T | 1 X 15 ml |
| MAL 25T | 1 X 7.5 ml |

PRINCIPLE OF THE METHOD

Microalbumin-turbilatex is a quantitative turbidimetric test for the measurement of microalbumin (μ ALB) in human urine. Latex particles coated with specific antibodies anti-human albumin are agglutinated when mixed with samples containing μ ALB. The agglutination causes an absorbance change, dependent upon the μ ALB contents of the patient sample that can be quantified by comparison from a calibrator of known μ ALB concentration.

CLINICAL SIGNIFICANCE

Microalbuminuria is at present defined as an excretion rate for albumin between 20 and 200 mg/L, which is already above normal values but still below the values seen in patients with "conventional" proteinuria.

Microalbuminuria is a marker of an increased risk of diabetic nephropathy as well as cardiovascular disease in patients with insulin-dependent diabetes mellitus as well as with non-insulin-dependent diabetes mellitus. More recently, microalbuminuria has been found to be associated with cardiovascular disease also in the non-diabetic population. In fact, microalbuminuria may show to be a risk factor of cardiovascular disease among otherwise apparently healthy people.

KIT CONTENT

| Reagent | 25 tests/kit | 50 tests/kit | Major ingredients |
|------------------------------|--------------|--------------|---|
| Reaction buffer (R1 reagent) | 1 x 6ml | 1 x 12ml | Glycine Buffer |
| R2 reagent | 1 x 1.5ml | 1 x 3 ml | Latex particles coated with goat IgG anti-human Albumin |
| Cuvette & Bead | 25 no's | 50 no's | / |
| Specification | 1 copy | 1 copy | / |
| Magnetic card | 1 piece | 1 piece | / |

SAMPLES COLLECTION AND PRESERVATION

24 hours or random/ first morning urine specimen. It is recommended to adjust the pH at 7.0 with NaOH/HCL 1 mol/L. Stable 7 days at 2-8°C when sodium azide 1 g/L is added to prevent contamination. Urine should be centrifuged before testing.

PRECAUTIONS

Components from human origin have been tested and found to be negative for the presence of HBsAg, HCV, and antibody to HIV (1/2). However handle cautiously as potentially infectious.

STORAGE AND STABILITY

All the components of the kit are stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use. Do not use reagents over the expiration date.

ASSAY PROCEDURE

1. THE REAGENT DOSE DISPENSING IN THE CUVETTE

| Reagent | Dose |
|---|-------------|
| Reaction Buffer (R1 Reagent) | 240 μ l |
| Sample | 5 μ l |
| anti-human Albumin - latex combo (R2 Reagent) | 60 μ l |

2. OPERATING STEP

- (1). When starting it shows "Read card", Put the corresponding lot magnetic card in the reader slot, reading the card correctly, the screening displays the reagent name and lot number, the instrument status indication light is working (yellow-green). Please check carefully.
- (2). After confirm the lot number, dispense 240 μ l reaction buffer (R1) to a colorimetric cuvette and 5 μ l sample and mixing sticker. Do not produce bubbles when dispensing the sample.
- (3). Put the colorimetric cuvette in the detection well, gently press the cuvette until it contacts the bottom. The status indication light will be off when the analyzer detects the colorimetric cuvette successfully
- (4). When the analyzer shows "Add [R2]", dispense 60 μ l anti-human Albumin antibody latex combo (R2) in the cuvette. The analyzer will mix automatically and start to detect, it shows: "Testing..." it shows the result automatically and record the value
- (5). After detecting, move the cuvette. The status indication light will be working (yellow- green). Return to step (2) to detect next sample.
- (6). If moving the cuvette out in the process of detecting, the screen shows "Give up testing".
- (7). If the result shows >230.00mg/L, may use 0.01 mol/L PBS to dilute the sample 1:5 (add 400 μ l PBS in 100 μ l sample), input dilution multiply 5, the analyser can calculate the sample concentration automatically.

QUALITY CONTROL

To ensure adequate quality control, Normal and abnormal control with assayed values should be run as unknown samples.

REFERENCE VALUES

Normal values up to 30 mg/24 hrs urine specimen and 20 mg/L in a first morning urine specimen.
Each laboratory should establish its own reference range.

INTERFERENCES

Glucose (2 g/L), hemoglobine (10 g/L) and creatinine (3 g/L), do not interfere. Urea (1 g/L) and bilirubin (10 mg/dL), interfere. Other substances may interfere.

BIBLIOGRAPHY

1. Feldt-Rasmussen B et al. J Diab Comp 1994; 8: 137-145.
2. Panuyiotou B N. Journal International Medical Research 1994; 22: 181-201.
3. Bar J et al. Diabetic Medicine 1995; 12: 649-656.
4. Gilbert R E et al. Diabetic Medicine 1994; 11: 636-645.
5. Medcalf E A et al. Clin Chem 1990; 36/3: 446-449.
6. Young DS. Effects of drugs on clinical laboratory test, 4th ed. AACC Press, 1995.